

**Date:****Initials:**

This document contains supply checklists and the protocol for the QuantiFluor dsDNA System using the GloMax Fluorometer to measure DNA quantity (concentration) prior to sample submission for Genotyping-by-Sequencing (GBS).

\* All solutions, standards and unknowns should be at room temperature prior to pipetting.

**Supplies Checklist:**

- Gloves
- Pipettes (list sizes)
  - Multi-channel
  -
- Pipette tips
- 96-well plate
- 1.5 ml microcentrifuge tubes
- Falcon tubes
- Dish thingy - does this have a specific name?
- USB-drive!

**Reagents Checklist:**

- Promega QuantiFluor dsDNA System kit (small yellow box in freezer)
  - 20X TE Buffer
  - QuantiFluor dsDNA dye
  - Lambda dsDNA
- Nuclease-free water
- Experimental ("unknown") DNA samples

**1. Prepare 1X TE Buffer:**

- Dilute 20X TE buffer 20-fold (1:20) with nuclease-free water. When doing a full plate, prepare enough diluted buffer for both serial dilutions and to create the working solution for 100 samples.
  - \*\* Remember: 1, 000 ul = 1 ml

<b>20X TE buffer</b>	<b>1000 ul + 30 ul</b>
<b><u>Nuclease-free water</u></b>	<b><u>19, 000 ul + 570 ul</u></b>
<b>1X TE buffer</b>	<b>20, 000 ul + 600 ul</b>

**2. Prepare Working Solution:**

- Dilute the QuantiFluor dsDNA dye 1:400 in 1X TE buffer to make the working solution. For a full plate, you need enough working solution for 200 ul in each of the 96 wells (19, 200 ul).
  - The working solution is stable for at least 2 hours at 25 C.

<b>QuantiFluor dsDNA dye</b>	<b>50 ul</b>
<b><u>1X TE buffer</u></b>	<b><u>19, 950 ul</u></b>
<b>Working solution</b>	<b>20, 000 ul</b>

**3. Prepare a Standard Curve:**

- Label seven in 1.5 ml tubes for each standard as follows:
  - 200, 50, 12.5, 3.1, 0.78, 0.2, and 0.05.
- Prepare dsDNA standards by mixing the Lambda DNA Standard (100 ng/ul) and serially diluting as shown in table below.
  - Shake DNA to mix since the concentration gradient. Do not vortex since DNA is fragile.
  - Take care not to introduce air bubbles. If this happens, spin them down.
  - Vortex well after each subsample - common error affecting standard curve.

<b>Preparing recommended dsDNA standard curve samples:</b>			
<b>Standard</b>	<b>Volume of dsDNA Standard</b>	<b>Volume of 1X TE Buffer (ul)</b>	<b>Final dsDNA Concentration (ng/ul)</b>
A	20 ul	80	20
B	25 ul of Standard A	75	5.0
C	25 ul of Standard B	75	1.25
D	25 ul of Standard C	75	0.31
E	25 ul of Standard D	75	0.078
F	25 ul of Standard E	75	0.02
G	25 ul of Standard F	75	0.005



